



TOTAL SYSTEM
ACCOUNTABILITY


Veterinary DR System



ClearVet System X

ClearVet Flat Panel Detector


Available in 17"x17" or 14"x17"
tethered or wireless



888.998.2838

info@clearvetdigital.com

www.clearvetdigital.com



4-way Float HF Table

Scratch resistant and contamination resistant table top for easy maintenance.

High Frequency Generator

42kW, 400mA, 125kVp

ClearVet FP

Image Field Sizes	17"x17" or 14"x17"
Spatial Resolution	3.9 lp/mm
Scintillator	CSI:TCWC
Number of Bits	14 bit

ClearVet System X Includes

- 4-way Float HF Table
- Cesium Flat Panel
- Acquisition Software, PC and monitor



Image Available in Seconds

Review patient images immediately. Adjust image and diagnose within seconds.

Global Support Team

With hundreds of digital systems installed in North America, ClearVet has extensive experience in veterinary radiography. Our team is here to support you.

Authorized Distributor



800.211.0124 • www.apexx360.com

For more information call 888.998.2838