



Tough Enough.

We passed all the tests. No one else has.

Engineered to perform in adverse and demanding conditions.
The ClearVet™CR* is the US Military's system of choice.

Independently Tested and Certified for Environmental Durability and Ruggedness

Some of the tests the ClearVet™CR* passed and others failed:

- VIBRATION TESTING—10 minutes per axis. PASSED
- SHOCK TESTING—30G Shock—18 Repetitions. PASSED
- TRANSIT DROP TEST—10 drops from two feet. PASSED
- BLOWING SAND TEST—Blowing Sand at 45 MPH for hours. PASSED
- SETTLING DUST—6 gm/m²/day for three days. PASSED
- HUMIDITY TEST—5 day extreme humidity cycle. PASSED
- SALT FOG TEST—5% Sodium Chloride Fog for 48 hours. PASSED
- FUNGUS TEST—28 day exposure to Fungus Spores. PASSED
- TEMPERATURE TESTING—Temperatures from of 30 degrees F to 105 degrees F. PASSED

The system remained operational through ALL tests

**Durability makes the ClearVet™CR the superior system for your
Mobile Veterinary application**

